Probe cover is a medical device used by healthcare provider. It is intended for barrier protection during invasive medical examination and diagnosis using transducer probe. Probe cover use to prevent cross-contamination of Transducer probe during endocavity (i.e. transvaginal or transrectal) ultrasound procedure.

For use, you should be trained in ultrasonography.

Intended user: It is intended to be used by healthcare provider trained in the use of ultrasound. The user should have no unique skills or user abilities required for protective probe cover use.

Target user group may include, but not limited to: Anesthesiologists, Department Nursing staff, Emergency medicine physicians, Endocrinologist, Interventional Radiologist & Interventional Radiology Technologist, Radiologist, Sonographers, Surgeons, Surgical Technicians and other clinicians.

Indication: For use as single use components as infectious waste.

Do NOT REUSE probe covers. Reuse or reprocessing may create a risk of contamination of the device and may cause patient infection or cross-infection.


Key BENEFITS : Prevent cross-contamination of probes and consequently reduce the risk of cross-contamination. For use in patients undergo invasive ultrasound examinations and diagnoses (i.e. transvaginal or transrectal scan)

Use ONLY WATER BASED LUBRICANT, avoid petroleum or mineral oil based lubricant that may harm the probe covers.

Disinfecting the transducer probe after each examination could reduce the probability of cross-contamination.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority.

This Product Contains Natural Rubber Latex Which May Cause Allergic Reactions.

Avoid damage to the cover by sharp edges, fingernails and etc. Place the cover over the transducer. Pull cover tightly over the transducer to remove wrinkles and air bubbles, taking care to avoid puncturing the cover. Inspect the cover to ensure no holes or tears.

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